

Acetyl-Histone H3 (K23) Rabbit Monoclonal Antibody (Q16695 / P68431) For Wb, If/Icc, Elisa

Item Number: CM0000184



Introduction

This rabbit monoclonal antibody specifically detects acetylated histone H3 at lysine 23 (H3K23ac), with wide cross-reactivity to human, mouse, and rat. Suitable for Western blot, IF/ICC, and ELISA applications in epigenetic research and chromatin dynamics studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K23 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; H3C1; Acetyl-Histone H3-K23
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic Peptide A synthetic acetylated peptide around K23 of human Histone H3 (P68431).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	16kDa
Shipping	Ice bag