

## Acetyl-Histone H3-K23 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Chip, Ihc-P, If/Icc, Elisa

Item Number: CM0000120



### Introduction

Acetyl-Histone H3-K23 rabbit polyclonal antibody (SWISS: Q16695/P68431) for Western blot, IHC-P, IF/ICC, ELISA, and ChIP. Cross-reacts with human, mouse, rat, and predicted wide-range species. Core nucleosome component with critical roles in chromatin regulation.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | Acetyl-Histone H3-K23 Rabbit pAb                                                                     |
| Remarks/Alias    | H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K23 |
| Species          | Human                                                                                                |
| GeneID (Human)   | H3                                                                                                   |
| GeneID           | 8290/8350                                                                                            |
| Immunogen        | A synthetic acetylated peptide around K23 of human Histone H3 (NP_003520.1)                          |
| Source           | Rabbit                                                                                               |
| Category         | Acetylated Antibodies                                                                                |
| Application      | WB, IHC-P, IF/ICC, ELISA, ChIP                                                                       |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                                      |
| SWISS            | Q16695/P68431                                                                                        |
| Protein Weight   | 15kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |