


Acetyl-Histone H3-K27 Rabbit Monoclonal Antibody For Chip, Wb & Ihc (Swiss: Q16695 / P68431)

Item Number: CM0000176



Introduction

Site-specific rabbit monoclonal antibody recognizing acetylated lysine 27 (K27ac) on histone H3 (SWISS Q16695/P68431). Validated for WB, ChIP, IHC-P, IF/ICC, ELISA, ChIP-seq. Wide predicted cross-reactivity (Human, Mouse, Rat). Suitable for epigenetic and chromatin biology research.

[Learn More](#)

Parameter	Details
Product Name	Acetyl-Histone H3-K27 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K27
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic acetylated peptide around K27 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag