

Acetyl-Histone H3-K27 Rabbit Monoclonal Antibody (Q16695 / P68431) For Wb, Ihc-P, Chip

Item Number: CM0000150



Introduction

Rabbit monoclonal antibody against acetylated lysine 27 of Histone H3 (Q16695/P68431). Suitable for WB, DB, IHC-P, IF/ICC, IP, ELISA, ChIP. Reacts with human, mouse, rat, and other species. Ideal for epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K27 Rabbit mAb
Species	Human
GeneID	8290/8350
Immunogen	Synthetic Peptide: A synthetic acetylated peptide around K27 of human Histone H3 (P68431)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag