


Acetyl-Histone H3-K27 Rabbit Monoclonal Antibody (Q16695/P68431) For Wb, Chip, If, Ihc

Item Number: CM0000197



Introduction

Acetyl-Histone H3-K27 Rabbit mAb (SWISS Q16695/P68431) validated for WB, ChIP, IHC-P, IF/ICC, ELISA, ChIP-seq, CUT&Tag. Cross-reacts with Human, Mouse, Rat. Ideal for histone modification and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K27 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K27
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic acetylated peptide around K27 of human Acetyl-Histone H3-K27 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq, CUT&Tag
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag