


Acetyl-Histone H3-K27 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Chip-Seq, Elisa

Item Number: CM0000200



Introduction

High-quality rabbit polyclonal antibody targeting acetylated Histone H3 at Lys27 (SWISS Q16695/P68431). Validated for multiple applications including WB, IHC-P, IF/ICC, IP, ELISA, ChIP, and ChIP-seq. Predicted broad species reactivity. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K27 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K27
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	A synthetic acetylated peptide around K27 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag