


Acetyl-Histone H3-K27 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Db, Elisa, Chip - Camelbio

Item Number: CM0000118



Introduction

High-quality rabbit polyclonal antibody specific to acetylated Histone H3 Lys27 (H3K27ac). Validated for WB, DB, ELISA, ChIP. Cross-reacts with human, mouse, rat. Ideal for epigenetics research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K27 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K27
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic acetylated peptide around K27 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag