

Acetyl-Histone H3-K36 (Q16695/P68431) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa

Item Number: CM0000146



Introduction

Rabbit polyclonal antibody targeting acetylated lysine 36 of human Histone H3 (Q16695/P68431). Validated for WB, IF/ICC, ELISA; predicted cross-reactivity with human, mouse, rat, and other species. Suitable for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K36 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K36
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic Peptide: A synthetic acetylated peptide around K36 of human Histone H3 (NP_003484.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag