

Acetyl-Histone H3 (K36) Rabbit Polyclonal Antibody (Q16695/P68431) For Wb And Elisa

Item Number: CM0000178



Introduction

Rabbit polyclonal antibody specific for acetylated lysine 36 of histone H3 (SWISS: Q16695/P68431). Suitable for Western blot and ELISA, with predicted wide species cross-reactivity. For epigenetic and chromatin research applications.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K36 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K36
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic acetylated peptide around K36 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag