


## Acetyl-Histone H3-K4 (Q16695/P68431) Rabbit Polyclonal Antibody For Wb And Elisa

Item Number: CM0000141



### Introduction

Rabbit polyclonal antibody specific for acetylated histone H3 at lysine 4 (K4). Suitable for Western blot and ELISA in human, mouse, and rat samples, with predicted wide range cross-reactivity. Ideal for epigenetic research and IVD development.

[Learn More](#)

| Parameter        | Value                                                                      |
|------------------|----------------------------------------------------------------------------|
| Product Name     | Acetyl-Histone H3-K4 Rabbit pAb                                            |
| Remarks/Alias    | H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K4                |
| Species          | Human                                                                      |
| GeneID (human)   | H3                                                                         |
| GeneID           | 8290/8350                                                                  |
| Immunogen        | A synthetic acetylated peptide around K4 of human Histone H3 (NP_003484.1) |
| Source           | Rabbit                                                                     |
| Category         | Acetylated Antibodies                                                      |
| Application      | WB, ELISA                                                                  |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                            |
| SWISS            | Q16695/P68431                                                              |
| Protein Weight   | 15 kDa                                                                     |
| Shipping         | Ice bag                                                                    |