


Acetyl-Histone H3-K4 Rabbit Monoclonal Antibody (Q16695/P68431) For Wb, Ip, Ihc, Chip, Elisa

Item Number: CM0000158



Introduction

Rabbit monoclonal antibody against acetylated Histone H3 at lysine 4 (Acetyl-Histone H3-K4). Validated for WB, IP, IF/ICC, IHC-P, DB, ChIP, ELISA. Cross-reacts with human, mouse, rat. Target weight ~15 kDa. Suitable for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K4 Rabbit mAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K4
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K4 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IP, IF/ICC, IHC-P, DB, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag