


Acetyl-Histone H3-K4 Rabbit Pab (Q16695 / P68431) For Wb, If, Elisa, Chip

Item Number: CM0000161



Introduction

Acetyl-Histone H3-K4 Rabbit Polyclonal Antibody (Q16695/P68431) detects endogenous acetyl-H3K4. Validated for WB, IF/ICC, ELISA, and ChIP in human, mouse, and rat, with predicted wide cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K4 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K4
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic Peptide - A synthetic acetylated peptide around K4 of human Histone H3 (NP_003484.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag