


Acetyl-Histone H3 (K4/K9/K14/K18/K23/K27) Rabbit Polyclonal Antibody (Q16695 / P68431) For Wb, Elisa

Item Number: CM0000166



Introduction

Polyclonal rabbit anti-acetyl-histone H3 antibody targeting K4, K9, K14, K18, K23, K27 modifications. Validated for WB and ELISA with broad cross-reactivity (human, mouse, rat). Ideal for epigenetics research and chromatin studies. 15 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K4/K9/K14/K18/K23/K27 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K4/K9/K14/K18/K23/K27
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic acetylated peptide around K4 & K27 of human Acetyl-Histone H3-K4/K9/K14/K18/K23/K27 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag