


Acetyl-Histone H3-K56 (Q16695 / P68431) Rabbit Monoclonal Antibody For Chip, Wb, If

Item Number: CM0000105



Introduction

Acetyl-Histone H3-K56 Rabbit mAb (clone A22565) recognizes acetylated lysine 56 of human Histone H3 (Q16695/P68431). Validated in ChIP-seq, WB, IF, and CUT&Tag for epigenetics research. Cross-reacts with mouse, rat, and a wide range of predicted species.

[Learn More](#)

Parameter	Specification
Product Name	Acetyl-Histone H3-K56 Rabbit mAb
Remarks/Alias	H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; Acetyl-Histone H3-K56
Species	Human
GeneID (human)	8290/ 8350
GeneID	8290/8350
Immunogen	Synthetic peptide; A synthetic acetylated peptide around K56 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IP, IF/ICC, IHC-P, DB, ChIP, ChIP-seq, CUT&Tag, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag