


Acetyl-Histone H3-K56 Rabbit Polyclonal Antibody (Q16695/P68431) For Epigenetics Research | Wb, Chip, Elisa

Item Number: CM0000203



Introduction

Acetyl-specific rabbit polyclonal antibody targeting histone H3 acetylated at lysine 56. Validated for Western blot, IF/ICC, ELISA, and ChIP, with broad cross-reactivity. Ideal for chromatin and epigenetic modification studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K56 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K56
Species	Human
GeneID (human)	H3C family genes; H3-4 (HIST3H3)
GeneID	8290/8350
Immunogen	A synthetic acetylated peptide around K56 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag