


Acetyl-Histone H3 K9 (Q16695/P68431) Rabbit Polyclonal Antibody - Wb Ip If Ihc Chip Elisa

Item Number: CM0000170



Introduction

Acetyl-Histone H3-K9 rabbit polyclonal antibody for detection in human, mouse, rat. Suitable for WB, IP, IF/ICC, IHC-P, ChIP, ChIP-seq, ELISA. 15 kDa, broad cross-reactivity.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K9 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K9
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IP, IF/ICC, IHC-P, ChIP, ChIP-seq, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag