

Acetyl-Histone H3-K9 Rabbit Mab (Q16695/P68431) For Epigenetics Research And Chip Applications

Item Number: CM0000172



Introduction

Rabbit monoclonal antibody against acetylated histone H3 lysine 9 (H3K9ac) for human, mouse, rat. Validated for WB, ChIP-seq, IF/ICC, and more. Ideal for epigenetic, chromatin, and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K9 Rabbit mAb
Remarks/Alias	H3A, HIST1H3A, H3C1, H3C2, H3C3, H3C4, H3C6, H3C7, H3C8, H3C10, H3C11, H3C12, Acetyl-Histone H3-K9
Species	Human
GeneID (Human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic acetylated peptide around K9 of human Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, IP, ELISA, ChIP, ChIP-seq, CUT&Tag
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag