

Acetyl-Histone H3 (K9/K14/K18/K23/K27) Rabbit Polyclonal Antibody For Wb, Elisa - Q16695/P68431

Item Number: CM0000149



Introduction

Rabbit polyclonal antibody targeting acetylated histone H3 at lysines 9, 14, 18, 23, and 27. Validated for WB and ELISA, with predicted wide cross-reactivity. Suitable for epigenetic research and chromatin biology studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K9/K14/K18/K23/K27 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K9/K14/K18/K23/K27
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic acetylated peptide around K9, K14, K18, K23, and K27 of human acetyl-Histone H3 (NP_003520.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag