

Acetyl-Histone H3-K9/K14/K18/K23/K27 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, If/Icc, Elisa

Item Number: CM0000135



Introduction

Rabbit polyclonal antibody targeting acetylated Histone H3 at lysines 9, 14, 18, 23, and 27. Validated for Western blot, immunofluorescence/ICC, and ELISA. Predicted cross-reactivity with human, mouse, rat, and other species. Suitable for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K9/K14/K18/K23/K27 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K9/K14/K18/K23/K27
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag