

Acetyl-Histone H3 (Lys4) Rabbit Monoclonal Antibody For Wb, Chip, If | Q16695/P68431

Item Number: CM0000147



Introduction

Rabbit monoclonal antibody specific for acetylated Histone H3 at Lys4 (Q16695/P68431). Suitable for Western blot, ChIP, IF/ICC, IP, DB, and ELISA. Predicted to cross-react with human, mouse, and rat samples. Ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H3-K4 Rabbit mAb
Remarks/Alias	H3.4; H3/g; H3FT; H3t; HIST3H3; Histone H3; HIST1H3A; Acetyl-Histone H3-K4
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide A synthetic acetylated peptide around K4 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IF/ICC, IP, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag