


Acetyl-Histone H4-K12 Rabbit Mab (P62805) For Chip, Wb, And Ihc-P | Camelbio

Item Number: CM0000195



Introduction

Rabbit monoclonal antibody specific for acetylated lysine 12 of human histone H4. Suitable for WB, IP, IHC-P, ChIP, and ELISA. Cross-reacts with human, mouse, rat, and other species. Ideal for epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K12 Rabbit mAb
Remarks/Alias	H4/I; H4C1; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FI; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4B; Acetyl-Histone H4-K12
Species	Human
GeneID (human)	8359/8370
GeneID	8359/8370
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K12 of human Histone H4 (NP_003486.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IP, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11 kDa
Shipping	Ice bag