


Acetyl-Histone H4 (K12) Rabbit Monoclonal Antibody For Wb, Chip, Chip-Seq, Elisa | P62805

Item Number: CM0000164



Introduction

Specific rabbit mAb targeting acetylated lysine 12 of human histone H4. Suitable for WB, DB, IHC-P, ELISA, ChIP, ChIP-seq. Cross-reacts with human, mouse, rat, and predicted wide-range species. Ideal for epigenetic and chromatin biology studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K12 Rabbit mAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Acetyl-Histone H4-K12
Species	Human
GeneID (human)	8359/8370
GeneID	8359
Immunogen	Synthetic acetylated peptide around K12 of human Histone H4 (P62805)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag