

Acetyl-Histone H4 (K12) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62805

Item Number: CM0000125



Introduction

Specific rabbit polyclonal antibody targeting acetylated lysine 12 of human Histone H4 (P62805). Validated for WB, IHC-P, IF/ICC, and ELISA with broad cross-reactivity to mouse, rat, and other species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K12 Rabbit pAb
Remarks/Alias	H4/I; H4C1; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FI; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4B; Acetyl-Histone H4-K12
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic peptide - A synthetic acetylated peptide around K12 of human Histone H4 (NP_001029249.1).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag