


Acetyl-Histone H4-K16 (P62805) Rabbit Polyclonal Antibody For Wb And Elisa - Epigenetics Research

Item Number: CM0000189



Introduction

Acetyl-Histone H4-K16 rabbit polyclonal antibody (UniProt P62805) validated for WB and ELISA. Broad cross-reactivity with human, mouse, and rat. Ideal for studying histone acetylation, chromatin dynamics, and epigenetic regulation in diagnostic research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K16 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Acetyl-Histone H4-K16
Species	Human
GeneID (human)	H4
GeneID	8370
Immunogen	A synthetic acetylated peptide around K16 of human Histone H4 (NP_003539.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag