


Acetyl-Histone H4-K16 Rabbit Monoclonal Antibody (P62805) For Wb, Chip, And Ihc-P Applications

Item Number: CM0000123



Introduction

Highly specific rabbit monoclonal antibody targeting acetylated lysine 16 of Histone H4 (SWISS P62805). Validated for WB, DB, IHC-P, ELISA, and ChIP with broad reactivity (human, mouse, rat). Essential for chromatin and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K16 Rabbit mAb
Remarks/Alias	H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FA; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4A; Acetyl-Histone H4-K16
Species	Human
GeneID (human)	8359/8370
GeneID	8359/8370
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K16 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11 kDa
Shipping	Ice bag