


Acetyl-Histone H4-K5 Rabbit Mab (P62805) For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000112



Introduction

Acetyl-Histone H4-K5 Rabbit mAb (P62805) detects acetylated Histone H4 at Lys5 in human, mouse, rat. Validated for WB, DB, IHC-P, IF/ICC, ELISA. Ideal for epigenetic research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K5 Rabbit mAb
Remarks/Alias	H4/p; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; HIST4H4; Acetyl-Histone H4-K5
Species	Human
GeneID (human)	8359
GeneID	8359
Immunogen	Synthetic Peptide: A synthetic acetylated peptide around K5 of human Histone H4 (P62805)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag