


Acetyl-Histone H4-K5 Rabbit Polyclonal Antibody (P62805) - Validated For Wb, Ihc-P, Elisa, Chip

Item Number: CM0000153



Introduction

Rabbit polyclonal antibody targeting acetylated Histone H4 at Lys5 (SWISS P62805). Suitable for WB, IHC-P, ELISA, and ChIP applications. Reacts with human, mouse, rat, and other species (wide range predicted). Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K5 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Acetyl-Histone H4-K5
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H4 (NP_001029249.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag