

Acetyl-Histone H4-K5 Rabbit Polyclonal Antibody (P62805) - Wb, Elisa, Chip Applications

Item Number: CM0000122



Introduction

Rabbit polyclonal antibody specific for acetylated lysine 5 of human Histone H4 (P62805). Validated for WB, ELISA, and ChIP, with wide predicted cross-reactivity among human, mouse, and rat. An essential tool for epigenetics and chromatin research.

[Learn More](#)

Parameters	Description
Product Name	Acetyl-Histone H4-K5 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Acetyl-Histone H4-K5
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic acetylated peptide around K5 of human Histone H4 (NP_001029249.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag