

Acetyl-Histone H4 (K5/K8/K12/K16) Rabbit Polyclonal Antibody P62805 For Wb, Ihc-P, Elisa

Item Number: CM0000133



Introduction

Rabbit polyclonal antibody recognizing acetylated Histone H4 at lysines 5, 8, 12, and 16 (SWISS P62805). Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat and other species. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K5/K8/K12/K16 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Acetyl-Histone H4-K5/K8/K12/K16
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	A synthetic acetylated peptide around K5/K8/K12/K16 of human Histone H4
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11 kDa
Shipping	Ice bag