


Acetyl-Histone H4 (K8) Rabbit Polyclonal Antibody (P62805) - Wb, Ihc-P, If/Icc, Elisa, Chip, Chip-Seq

Item Number: CM0000138



Introduction

Acetyl-Histone H4 (Lys8) rabbit polyclonal antibody (P62805) detects endogenous acetyl-histone H4 in human, mouse, rat, and other species. Validated for WB, IHC-P, IF/ICC, ELISA, ChIP, and ChIP-seq applications.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K8 Rabbit pAb
Remarks/Alias	H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FA; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4A; Acetyl-Histone H4-K8
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic peptide A synthetic acetylated peptide around K8 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag