


Acetyl-Histone H4 (K8Ac) Rabbit Monoclonal Antibody - P62805 | Validated In Wb, Ip, Ihc-P, Chip

Item Number: CM0000140



Introduction

Acetyl-Histone H4-K8 rabbit monoclonal antibody (P62805) detects acetylated lysine 8 of histone H4 in human, mouse, rat. Suitable for WB, IP, IHC-P, DB, ChIP, ELISA. Ideal for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K8 Rabbit mAb
Remarks/Alias	H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4FA; H4-16; H4C11; H4C12; H4C13; H4C14; H4C15; H4C16; HIST1H4A
Species	Human
GeneID (Human)	8359
Immunogen	Synthetic peptide - A synthetic acetylated peptide around K8 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IP, IHC-P, DB, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11 kDa
Shipping	Ice bag