

Acetyl-Histone H4 K91 Rabbit Polyclonal Antibody (P62805) - Wb, If/Icc, Elisa

Item Number: CM0000143



Introduction

Rabbit polyclonal antibody targeting acetylated lysine 91 of human Histone H4 (P62805). Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat, and other species. Ideal for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-Histone H4-K91 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Acetyl-Histone H4-K91
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic peptide: A synthetic acetylated peptide around K91 of human Histone H4 (NP_003529.1)
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag