


Acetyl-Hmgb1 (Lys29) Rabbit Polyclonal Antibody (P09429) For Wb, If/Icc, Elisa

Item Number: CM0000142



Introduction

High-quality rabbit polyclonal antibody specific for acetylated HMGB1 at Lys29 (SWISS P09429). Validated for WB, IF/ICC, ELISA with human, mouse, rat reactivity. Ideal for inflammation, DAMP signaling, and post-translational modification research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-HMGB1-K29 Rabbit pAb
Remarks/Alias	HMG1; HMG3; HMG-1; SBP-1; Acetyl-HMGB1-K29
Species	Human
GeneID (human)	3146
GeneID	3146
Immunogen	Synthetic Peptide A synthetic acetylated peptide around K29 of human HMGB1.
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09429
Protein Weight	25kDa
Shipping	Ice bag