

Acetyl-P53 (K370) Rabbit Monoclonal Antibody (P04637) For Western Blot And Elisa

Item Number: CM0000192



Introduction

Rabbit monoclonal antibody detecting acetylated human p53 at lysine 370 (SWISS P04637). Validated for Western blot and ELISA in human and mouse samples. Ideal for cancer and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-p53- K370 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Acetyl-p53- K370
Species	Human
GeneID (human)	7157
GeneID	7157
Immunogen	Synthetic Peptide A synthetic acetylated peptide around K370 of human p53 (P04637).
Source	Rabbit
Category	Acetylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag