

Acetyl-P53-K382 Rabbit Polyclonal Antibody (P04637) For Dot Blot And Elisa

Item Number: CM0000103



Introduction

Rabbit polyclonal antibody specific for acetylated lysine 382 of human p53 (P04637). Validated for dot blot and ELISA, it is ideal for studying p53 acetylation in cancer and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Acetyl-p53-K382 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Acetyl-p53-K382
Species	Human
GeneID (human)	7157
GeneID	7157
Immunogen	A synthetic acetylated peptide around K382 of human p53 (NP_000537.3)
Source	Rabbit
Category	Acetylated Antibodies
Application	DB, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag