

Acot9 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y305

Item Number: CM0013421

Introduction

ACOT9 Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA. Targets mitochondrial acyl-CoA thioesterase 9 in human, mouse, and rat samples. Ideal for studying lipid and amino acid metabolism.

[Learn More](#)

Parameter	Value
Product Name	ACOT9 Rabbit pAb
Remarks/Alias	ACATE2; CGI-16; MTACT48; MT-ACT48; ACOT9
Species	Human
Gene ID (Human)	23597
Gene ID	23597
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-250 of human ACOT9 (NP_001028755.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y305
Protein Weight	50kDa
Shipping	Ice bag