

Acsl4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc - O60488

Item Number: CM0021976

Introduction

High-quality ACSL4 Rabbit Polyclonal Antibody (O60488), validated for Western blot, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for lipid metabolism and ferroptosis research.

[Learn More](#)

Parameter	Value
Product Name	ACSL4 Rabbit pAb
Remarks/Alias	ACS4; FACL4; LACS4; MRX63; MRX68; XLID63; ACSL4
Species	Human
GeneID (Human)	2182
GeneID	2182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ACSL4 (NP_004449.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60488
Protein Weight	79kDa
Shipping	Ice bag