

Acsl5 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9ULc5

Item Number: CM0008808

Introduction

Highly specific ACSL5 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, and ELISA. Targets human long-chain-fatty-acid--CoA ligase 5 (LACS5), a key enzyme in fatty acid activation and metabolism. Ideal for research in lipid biology, cancer, and metabolic disorders.

[Learn More](#)

Parameter	Value
Product Name	ACSL5 Rabbit mAb
Remarks/Alias	ACS2, ACS5, FACL5
Species	Human
Gene ID (Human)	51703
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-683 of human ACSL5 (NP_976313.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9ULC5
Protein Weight	76kDa
Shipping	Ice bag