


Actc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P68032

Item Number: CM0026770



Introduction

High-quality ACTC1 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA in human, mouse, and rat samples. Target is cardiac alpha-actin, a key cytoskeletal protein. Ideal for cardiac muscle research.

[Learn More](#)

Parameter	Value
Product Name	ACTC1 Rabbit pAb
Remarks/Alias	ACTC; ASD5; CMD1R; CMH11; LVNC4; ACTC1
Species	Human
GeneID (Human)	70
GeneID	70
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ACTC1 (NP_005150.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68032
Protein Weight	42kDa
Shipping	Ice bag