


Activator Protein 2 (Ap-2/Tfap2A) Rabbit Pab - P05549

Item Number: CM0024188



Introduction

Anti-TFAP2A rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and rat samples. Targets transcription factor AP-2-alpha, involved in development and gene regulation. Ideal for IVD research.

[Learn More](#)

Parameter	Specification
Product Name	Activator protein 2 (AP-2/TFAP2A) Rabbit pAb
Remarks/Alias	AP-2; BOFS; AP2TF; TFAP2; AP-2alpha; Activator protein 2 (AP-2/TFAP2A)
Species	Human
GeneID (human)	7020
Immunogen	A synthetic peptide corresponding to amino acids 100-200 of human TFAP2A (NP_003211.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P05549
Protein Weight	48kDa
Shipping	Ice bag