

Actr1B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P42025

Item Number: CM0019686

Introduction

Rabbit polyclonal antibody targeting human beta-actin (ACTR1B). Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Suitable for studies of microtubule-based vesicle transport and centrosome function.

[Learn More](#)

Parameter	Value
Product Name	ACTR1B Rabbit pAb
Remarks/Alias	PC3; ARP1B; CTRN2; ACTR1B
Species	Human
Gene ID (Human)	10120
Gene ID	10120
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human ACTR1B (NP_005726.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42025
Protein Weight	42kDa
Shipping	Ice bag