


Acvr1 Rabbit Monoclonal Antibody For Wb And Elisa - Q04771

Item Number: CM0003970



Introduction

High-quality ACVR1 Rabbit monoclonal antibody validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat ACVR1 (activin receptor type-1, ALK2). Ideal for bone morphogenetic protein signaling pathway and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	ACVR1 Rabbit mAb
Remarks/Alias	FOP; ALK2; SKR1; TSRI; ACTRI; ACVR1A; ACVRLK2; ACVR1
Species	Human
GeneID	90
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 410-509 of human ACVR1 (Q04771).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04771
Protein Weight	57kDa
Shipping	Ice bag