

Adam19 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H013

Item Number: CM0014932

Introduction

Rabbit polyclonal antibody against ADAM19 (MLTNB, MADDAM), suitable for WB and ELISA. Cross-reacts with human, mouse, rat. Target is a transmembrane metalloprotease involved in neurogenesis and osteoblast differentiation.

[Learn More](#)

Parameter	Value
Product Name	ADAM19 Rabbit pAb
Remarks/Alias	MLTNB; FKSG34; MADDAM; ADAM19
Species	Human
GeneID (human)	8728
GeneID	8728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 203-460 of human ADAM19 (NP_150377.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H013
Protein Weight	105kDa
Shipping	Ice bag