

Adam9 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13443

Item Number: CM0027413

Introduction

Validate ADAM9 with this polyclonal antibody for Western blot, IF/ICC, and ELISA. Suitable for human, mouse, rat, and monkey samples. Ideal for tumorigenesis and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	ADAM9 Rabbit pAb
Remarks/Alias	MCMP; MDC9; CORD9; Mltng; ADAM9
Species	Human
GeneID (Human)	8754
GeneID	8754
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-210 of human ADAM9 (NP_003807.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	Q13443
Protein Weight	91 kDa
Shipping	Ice bag