

Adamts13 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q76Lx8

Item Number: CM0010865



Introduction

Recombinant ADAMTS13 rabbit monoclonal antibody suitable for Western blot, immunofluorescence, and ELISA. Recognizes human, mouse, and rat species. Useful for studying von Willebrand factor cleavage and thrombotic disorders.

[Learn More](#)

Parameter	Description
Product Name	ADAMTS13 Rabbit mAb
Remarks/Alias	VWF-CP; C9orf8; vWF-CP; ADAM-TS13; ADAMTS-13; ADAMTS13
Species	Human
Gene ID	11093
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human ADAMTS13 (Q76LX8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q76LX8
Protein Weight	154kDa
Shipping	Ice bag