


Adora2A Rabbit Pab - P29274

Item Number: CM0020246



Introduction

ADORA2A Rabbit pAb is a polyclonal antibody targeting human adenosine receptor A2a. Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat, supporting G protein-coupled receptor research.

[Learn More](#)

Parameter	Value
Product Name	ADORA2A Rabbit pAb
Remarks/Alias	A2aR; RDC8; ADORA2; ADORA2A
Species	Human
Gene ID (Human)	135
Gene ID	135
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ADORA2A (NP_000666.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29274
Protein Weight	45kDa
Shipping	Ice bag