

Adprm Rabbit Pab - Q3Lie5

Item Number: CM0014320

Introduction

ADPRM Rabbit pAb is a polyclonal antibody for Western blot (WB) and ELISA applications, targeting human ADPRM (MDS006/C17orf48) with cross-reactivity to mouse and rat. Ideal for studying ADP-ribose metabolism and immune signaling.

[Learn More](#)

Parameter	Value
Product Name	ADPRM Rabbit pAb
Remarks/Alias	MDS006; C17orf48; NBLA03831; ADPRM
Species	Human
Gene ID (Human)	56985
Gene ID	56985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 203-342 of human ADPRM (NP_064618.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q3LIE5
Protein Weight	40kDa
Shipping	Ice bag