

Afap1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N556

Item Number: CM0027030

Introduction

AFAP1 Rabbit Polyclonal Antibody validated for Western blot, immunofluorescence, and ELISA. Detects human, mouse, rat AFAP1 (81 kDa). Useful for studying actin dynamics, lamellipodia, and cancer cell migration.

[Learn More](#)

Parameter	Value
Product Name	AFAP1 Rabbit pAb
Remarks/Alias	AFAP; AFAP110; AFAP-110; AFAP1
Species	Human
Gene ID (Human)	60312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 461-730 of human AFAP1 (NP_940997.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N556
Protein Weight	81kDa
Shipping	Ice bag