

Affinity Gel-Conjugated Pan Succinyl-Lysine Rabbit Monoclonal Antibody For Ip-MS

Item Number: CM0000044

Introduction

CamelBio supplies affinity gel-conjugated Pan Succinyl-Lysine rabbit monoclonal antibody for IP-MS, enabling species-independent succinyl-lysine detection in post-translational modification studies. Ideal for raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	Affinity Gel-conjugated Pan Succinyl-Lysine Rabbit mAb
Immunogen	Synthetic peptide
Source	Rabbit
Category	Other Modified Antibodies
Application	IP-MS
Cross Reactivity	Species independent
Shipping	Ice bag