


Agfg1 Rabbit Pab - P52594

Item Number: CM0015275



Introduction

Rabbit polyclonal antibody targeting human AGFG1, validated for WB and ELISA. Cross-reacts with mouse. Useful for studying vesicle docking, RNA trafficking, and HIV-1 Rev cofactor function.

[Learn More](#)

Parameter	Value
Product Name	AGFG1 Rabbit pAb
Remarks/Alias	HRB; RAB; RIP; AGFG1
Species	Human
GeneID (Human)	3267
GeneID	3267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human AGFG1 (NP_004495.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52594
Protein Weight	58kDa
Shipping	Ice bag